


Anti-Human Cd70 (Cd27L) Monoclonal Antibody For Flow Cytometry - P32970

Item Number: CM0008162



Introduction

-conjugated rabbit monoclonal antibody targeting human CD70 (CD27L, TNFSF7), validated for flow cytometry. Suitable for immune cell interaction and CD27/CD70 pathway research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD70 mAb
Alias	CD27L; LPFS3; CD27-L; CD27LG; TNFSF7; TNLG8A
Species	Human
Gene ID	970
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 39-193 of human CD70 (P32970)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (Flow Cytometry)
Cross Reactivity	Human
Swiss-Prot	P32970
Protein Weight	21 kDa
Shipping	Ice bag