

Anti-Human Cd72 Monoclonal Antibody () For Flow Cytometry - P21854

Item Number: CM0012059

Introduction

High-quality -conjugated anti-human CD72 rabbit monoclonal antibody, specifically validated for flow cytometry applications. Targets B-cell differentiation antigen CD72, a co-receptor modulating B-cell receptor signaling and immune homeostasis.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD72 mAb
Remarks/Alias	CD72; CD72b; LYB2; CD72 molecule
Species	Human
Gene ID (Human)	971
Gene ID	971
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 117-359 of human CD72(NP_001773.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P21854
Protein Weight	40kDa
Shipping	Ice bag