

Anti-Human Cd79A Monoclonal Antibody For Flow Cytometry (Intracellular) - P11912

Item Number: CM0008356



Introduction

-conjugated rabbit monoclonal antibody targeting human CD79a. Suitable for intracellular flow cytometry to detect B-cell receptor complex alpha chain in human samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD79a mAb
Remarks/Alias	CD79A; IGA; MB-1; CD79a molecule
Species	Human
Gene ID (Human)	973
Gene ID	973
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 165-225 of human CD79a (NP_001774.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P11912
Protein Weight	20kDa/25kDa
Shipping	Ice bag