

Anti-Human Cd79A Monoclonal Antibody For Intracellular Flow Cytometry - P11912

Item Number: CM0004529



Introduction

Rabbit anti-human CD79a monoclonal antibody conjugated with for intracellular flow cytometry. Targets B-cell antigen receptor-associated protein alpha chain (MB-1, Ig-alpha). Validated in human samples. UniProt P11912.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD79a mAb
Remarks/Alias	IGA; MB1; MB-1; IGAAlpha
Species	Human
GeneID (Human)	973
GeneID	973
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 165-225 of human CD79a (NP_001774.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P11912
Protein Weight	25kDa
Shipping	Ice bag