

Anti-Human Cd90/Thy1 Rabbit Monoclonal Antibody For Flow Cytometry - P04216

Item Number: CM0008652

Introduction

-conjugated rabbit anti-human CD90/Thy1 monoclonal antibody for flow cytometry. Recombinant immunogen (aa 20-130). Validated for human samples. UniProt P04216.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD90/Thy1 mAb
Remarks/Alias	CD90; CDw90
Species	Human
Gene ID (Human)	7070
Gene ID	7070
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P04216
Protein Weight	18kDa
Shipping	Ice bag