

Anti-Human Cd94 Monoclonal Antibody For Fc - Q13241

Item Number: CM0011756



Introduction

-conjugated rabbit anti-human CD94 (KLRD1) monoclonal antibody for flow cytometry. Detects human CD94, an NK cell receptor involved in immune surveillance and tolerance. Suitable for immuno-oncology and infectious disease research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD94 mAb
Remarks/Alias	KLRD1; CD94; natural killer cells antigen CD94
Species	Human
GeneID (Human)	3824
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-179 of human CD94 (NP_002253.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q13241
Protein Weight	17kDa/20kDa
Shipping	Ice bag