

Anti-Human Cd96/Tactile Monoclonal Antibody For Flow Cytometry - P40200

Item Number: CM0005492



Introduction

High-quality rabbit monoclonal antibody against human CD96/TACTILE, validated for flow cytometry. Ideal for studying immune cell adhesion, T-cell and NK-cell activation in inflammation. - P40200

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD96/TACTILE mAb
Remarks/Alias	TACTILE
Species	Human
Gene ID (Human)	10225
Gene ID	10225
Immunogen	Recombinant protein corresponding to a sequence within amino acids 25-519 of human CD96/TACTILE (NP_937839.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P40200
Protein Weight	65kDa
Shipping	Ice bag