


Anti-Human Cd98 Monoclonal Antibody () For Flow Cytometry - P08195

Item Number: CM0012135



Introduction

ABflo® 594-conjugated rabbit monoclonal antibody against human CD98 (SLC3A2), verified for flow cytometry. Suited for immunology and cancer research, enabling expression profiling of amino acid transporter heavy chain.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD98 mAb
Remarks/Alias	4F2; CD98; MDU1; 4F2HC; 4T2HC; NACAE; CD98HC
Species	Human
GeneID (human)	6520
GeneID	6520
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 206-630 of human CD98 (NP_002385.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P08195
Protein Weight	57kDa/61kDa/67kDa/68kDa/71kDa
Shipping	Ice bag