

Anti-Human Cd99 Monoclonal Antibody () For If/Icc, Fc - P14209

Item Number: CM0008577



Introduction

Rabbit anti-human CD99 monoclonal antibody conjugated to . Suitable for immunofluorescence and flow cytometry. Validated for human samples. Ideal for T-cell adhesion and leukocyte extravasation research. SWISS: P14209.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD99 mAb
Remarks/Alias	MIC2; HBA71; MIC2X; MIC2Y; MSK5X
Species	Human
GeneID (Human)	4267
GeneID	4267
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 24-123 of human CD99 (NP_002405.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Human
SWISS	P14209
Protein Weight	19kDa
Shipping	Ice bag