

Anti-Human Cdcp1 Monoclonal Antibody For Flow Cytometry - Q9H5V8

Item Number: CM0011370

Introduction

-conjugated rabbit monoclonal antibody targeting human CDCP1 (CD318, TRASK, SIMA135), validated for flow cytometry. Used in cell adhesion, cancer metastasis, leukemia stem cell research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CDCP1 mAb
Remarks/Alias	CD318; TRASK; SIMA135
Species	Human
GeneID	64866
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-666 of human CDCP1 (NP_073753.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9H5V8
Protein Weight	37kDa/72kDa/92kDa
Shipping	Ice bag