


Anti-Human Ceacam6 Monoclonal Antibody (Conjugate) For Flow Cytometry - P40199

Item Number: CM0004046



Introduction

-conjugated rabbit monoclonal antibody targeting human CEACAM6 (CD66c). Validated for flow cytometry. Ideal for tumor biology, neutrophil adhesion, and cell migration studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CEACAM6 mAb
Remarks/Alias	NCA; CEAL; CD66c
Species	Human
GeneID	4680
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-320 of human CEACAM6 (NP_002474.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P40199
Protein Weight	37kDa
Shipping	Ice bag