

Anti-Human Clec12A/Cd371 Monoclonal Antibody For Flow Cytometry - Q5Qgz9-2

Item Number: CM0010303

Introduction

-conjugated rabbit monoclonal antibody targeting human CLEC12A/CD371 (CLL-1, MICL). Validated for flow cytometry and reacts with human samples. Suitable for immunophenotyping of myeloid cells.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CLEC12A/CD371 mAb
Remarks/Alias	CLL1; MICL; CD371; CLL-1; DCAL-2
Species	Human
GeneID	160364
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 65-265 of human CLEC12A/CD371(NP_612210.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q5QGZ9-2
Protein Weight	31KDa
Shipping	Ice bag