

Anti-Human Clec12A/Cd371 Monoclonal Antibody (Abflo 488) For Flow Cytometry - Q5Qgz9-2

Item Number: CM0004203

Introduction

ABflo 488-conjugated rabbit monoclonal antibody against human CLEC12A/CD371 (CLL-1, MICL). Designed for flow cytometry, targeting isoform 2 (aa 65-265). Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CLEC12A/CD371 mAb
Remarks/Alias	CLL1; MICL; CD371; CLL-1; DCAL-2
Species	Human
GeneID (Human)	160364
GeneID	160364
Immunogen	Recombinant fusion protein containing amino acids 65-265 of human CLEC12A/CD371 (NP_612210.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q5QGZ9-2
Protein Weight	8kDa/24kDa/27kDa/30kDa/32kDa
Shipping	Ice bag