

Anti-Human Cytokeratin 18 Monoclonal Antibody For Intracellular Flow Cytometry - P05783

Item Number: CM0004157

Introduction

-conjugated rabbit monoclonal antibody targeting human cytokeratin 18 (K18/CK-18/CYK18) for intracellular flow cytometry. Detects human K18 (48 kDa) in epithelial cells; suitable for cancer research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Cytokeratin 18 mAb
Remarks/Alias	K18; CK-18; CYK18
Species	Human
Gene ID (Human)	3875
Gene ID	3875
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 241-430 of human Cytokeratin 18 (NP_000215.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P05783
Protein Weight	48kDa
Shipping	Ice bag