


Abflo 594 Anti-Human Cytokeratin 19 (Krt19) Monoclonal Antibody For Intracellular Flow Cytometry - P08727

Item Number: CM0008954



Introduction

ABflo 594-conjugated rabbit monoclonal antibody targeting human Cytokeratin 19 (CK-19, KRT19). Designed for intracellular flow cytometry, reacting with human samples. Detects CK-19 at costameres and intermediate filaments. Protein weight 44kDa.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Cytokeratin 19 (KRT19) mAb
Remarks/Alias	K19; CK19; K1CS
Species	Human
Gene ID (Human)	3880
Gene ID	3880
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 301-400 of human Cytokeratin 19 (NP_002267.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P08727
Protein Weight	44kDa
Shipping	Ice bag