

## Anti-Human DLL4 Monoclonal Antibody ( ) For Fc - Q9Nr61

Item Number: CM0004043



### Introduction

-conjugated rabbit anti-human DLL4 monoclonal antibody for flow cytometry. Targets Notch ligand DLL4, involved in angiogenesis and neurogenesis. Validated for human samples. Ideal for IVD research and diagnostic assay development.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human DLL4 mAb                                                                                       |
| Remarks/Alias    | AOS6; delta4; hdelta2; delta-like protein 4                                                                      |
| Species          | Human                                                                                                            |
| GeneID (Human)   | 54567                                                                                                            |
| GeneID           | 54567                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 27-524 of human DLL4 (NP_061947.1) |
| Source           | Rabbit                                                                                                           |
| Category         | Monoclonal Antibodies                                                                                            |
| Application      | FC                                                                                                               |
| Cross Reactivity | Human                                                                                                            |
| SWISS            | Q9NR61                                                                                                           |
| Protein Weight   | 75kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |