

Anti-Human DLL4 Monoclonal Antibody () For Flow Cytometry - Q9Nr61

Item Number: CM0004065



Introduction

-conjugated rabbit monoclonal antibody targeting human DLL4 (Delta-like protein 4). Validated for flow cytometry detection of Notch ligand involved in angiogenesis and neurogenesis. Expressed in vascular endothelium.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human DLL4 mAb
Remarks/Alias	AOS6; delta4; hdelta2; delta-like protein 4
Species	Human
GeneID (Human)	54567
GeneID	54567
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-524 of human DLL4 (NP_061947.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9NR61
Protein Weight	75kDa
Shipping	Ice bag