

Anti-Human Eg-Vegf/Pk1 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P58294

Item Number: CM0003389

Introduction

High-quality rabbit monoclonal antibody targeting human EG-VEGF/PK1 (PROK1). Validated for WB, IF-P, ELISA; cross-reacts with mouse and rat. Ideal for angiogenesis and reproductive biology research.

[Learn More](#)

Parameter	Value
Product Name	Human EG-VEGF/PK1 Rabbit mAb
Remarks/Alias	PK1; PRK1; EGVEGF; PROK1
Species	Human
Gene ID (Human)	84432
Gene ID	84432
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-105 of human PROK1 (P58294)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P58294
Protein Weight	12kDa
Shipping	Ice bag