

Anti-Human Fap Monoclonal Antibody For Fc - Q12884

Item Number: CM0011741



Introduction

-conjugated rabbit monoclonal antibody against human FAP, validated for flow cytometry. FAP is a serine protease involved in extracellular matrix remodeling and tumor progression. Ideal for cancer and fibrosis research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human FAP mAb
Remarks/Alias	FAP; DPPIV; FAPA; FAPalpha; SIMP; prolyl endopeptidase FAP
Species	Human
GeneID (human)	2191
GeneID	2191
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q12884
Protein Weight	27kDa/87kDa
Shipping	Ice bag