

Anti-Human Fatp1 () Rabbit Monoclonal Antibody For Flow Cytometry - Q6Pcb7

Item Number: CM0004338

Introduction

-conjugated rabbit monoclonal antibody targeting human long-chain fatty acid transport protein 1 (FATP1/SLC27A1), suitable for intracellular flow cytometry (FC). Enables research into fatty acid import, lipid metabolism, and related pathways.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human FATP1 mAb
Remarks/Alias	FATP; FATP1; ACSVL5; FATP-1
Species	Human
GeneID (human)	376497
GeneID	376497
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 510-646 of human FATP1 (NP_940982.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	Q6PCB7
Protein Weight	52kDa/71kDa
Shipping	Ice bag