


Anti-Human FcεR1A Monoclonal Antibody () For Fc - P12319

Item Number: CM0005632



Introduction

Rabbit anti-Human FcεR1α monoclonal antibody, -conjugated, validated for flow cytometry detection of the high-affinity IgE receptor (FcεR1α) on eosinophils, mast cells, and basophils in allergy and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human FcεR1α mAb
Remarks/Alias	FCE1A, FcERI, FCERIA
Species	Human
GeneID (Human)	2205
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-205 of human FcεR1α (NP_001992.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P12319
Protein Weight	30kDa
Shipping	Ice bag