


Anti-Human FcεR1A Monoclonal Antibody For Fc - P12319

Item Number: CM0004254



Introduction

High affinity IgE receptor alpha (FcεR1α, FCER1A) rabbit monoclonal antibody conjugated to for flow cytometry (FC) detection in human samples. Recombinant immunogen from amino acids 26-205.

[Learn More](#)

Parameter	Details
Product Name	Rabbit anti-Human FcεR1α mAb
Remarks/Alias	FCE1A; FcERI; FCERIA
Species	Human
GeneID (Human)	2205
GeneID	2205
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P12319
Protein Weight	30kDa
Shipping	Ice bag