

Anti-Human Fscn2 (Fascin-2) Rabbit Polyclonal Antibody For Wb And Elisa - O14926

Item Number: CM0022040

Introduction

High-quality rabbit polyclonal antibody targeting human FSCN2 (Fascin-2), suitable for WB and ELISA. Cross-reactive with mouse and rat. Ideal for actin bundling, photoreceptor disk morphogenesis, and retinal degeneration studies.

[Learn More](#)

Parameter	Value
Product Name	FSCN2 Rabbit pAb
Remarks/Alias	RFSN; RP30; FSCN2
Species	Human
GeneID (Human)	25794
GeneID	25794
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-310 of human FSCN2 (NP_036550.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O14926
Protein Weight	55kDa
Shipping	Ice bag