


Anti-Human Galectin-3 Monoclonal Antibody For Flow Cytometry (Intracellular) - P17931

Item Number: CM0009263



Introduction

High-quality rabbit monoclonal antibody against human Galectin-3 for intracellular flow cytometry. Galectin-3 is a galactose-specific lectin involved in inflammation, splicing, and autophagy, making it a key diagnostic target.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Galectin 3 mAb
Remarks/Alias	L31; GAL3; MAC2; CBP35; GALBP; GALIG; LGALS2
Species	Human
GeneID (human)	3958
GeneID	3958
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human Galectin 3 (NP_002297.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P17931
Protein Weight	26kDa
Shipping	Ice bag