

Abflo 647 Anti-Human Gp6/Gpvi Monoclonal Antibody For Flow Cytometry - Q9Hcn6

Item Number: CM0011782



Introduction

ABflo 647-conjugated rabbit monoclonal antibody targeting human platelet glycoprotein VI (GP6/GPVI) for flow cytometry. Involved in collagen-induced platelet activation and thrombus formation. Validated for human samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human GP6/GPVI mAb
Remarks/Alias	GPV; GPVI; BDPLT11
Species	Human
GeneID (Human)	51206
GeneID	51206
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-267 of human GP6/GPVI (NP_057447.5)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9HCN6
Protein Weight	37kDa
Shipping	Ice bag