

Anti-Human Gp6/Gpvi Monoclonal Antibody For Flow Cytometry - Q9Hcn6

Item Number: CM0008945



Introduction

-conjugated rabbit monoclonal antibody against human GP6/GPVI, a key platelet collagen receptor. Validated for flow cytometry, reacts specifically with human GPVI. Suitable for thrombosis, hemostasis, and platelet activation studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human GP6/GPVI mAb
Remarks/Alias	GPV; GPVI; BDPLT11
Species	Human
Gene ID (Human)	51206
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9HCN6
Protein Weight	37kDa
Shipping	Ice bag