


Anti-Human Gprc5D Monoclonal Antibody For Flow Cytometry - Q9Nzd1

Item Number: CM0007724



Introduction

-conjugated rabbit monoclonal antibody targeting human GPCR5D, validated for flow cytometry. Detects a GPCR involved in hair and nail development. Ideal for immunodetection in human samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human GPCR5D mAb
Remarks/Alias	GPCR5D
Species	Human
Gene ID (Human)	55507
Gene ID	55507
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human GPCR5D (NP_061124.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9NZD1
Protein Weight	34kDa/36kDa/38kDa
Shipping	Ice bag