


Anti-Human Hbg1/2 Rabbit Monoclonal Antibody For Wb, If-P, And Elisa - P69891 / P69892

Item Number: CM0002509



Introduction

Rabbit monoclonal antibody targeting human HBG1/2 (gamma-1 and gamma-2 globin chains). Suitable for WB, IF-P, and ELISA. Ideal for fetal hemoglobin studies, hemoglobinopathy research, and IVD test development.

[Learn More](#)

Parameter	Value
Product Name	HBG1/2 Rabbit mAb
Remarks/Alias	HBG-T2; HBGA; HBGR; HSGGL1; PRO2979; HBG1
Species	Human
GeneID (Human)	3047/3048
GeneID	3047/3048
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 68-147 of human HBG1 (P69891).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human
SWISS	P69891/P69892
Protein Weight	16kDa
Shipping	Ice bag