

Anti-Human Iars (Isoleucyl-Trna Synthetase) Polyclonal Antibody For Wb, Elisa - P41252

Item Number: CM0015692

Introduction

Rabbit polyclonal antibody against human cytoplasmic isoleucyl-tRNA synthetase (IARS/IRS). Validated for Western blot and ELISA; detects a 144 kDa protein. Suitable for protein biosynthesis and translation studies.

[Learn More](#)

Parameter	Value
Product Name	IARS Rabbit pAb
Remarks/Alias	IRS; IARS; ILRS; ILERS; GRIDHH; PRO0785
Species	Human
Gene ID (Human)	3376
Gene ID	3376
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1158-1262 of human IARS (NP_002152.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P41252
Protein Weight	144kDa
Shipping	Ice bag