


Anti-Human Icam-1/Cd54 Monoclonal Antibody For Flow Cytometry - P05362

Item Number: CM0005296



Introduction

-conjugated rabbit anti-human ICAM-1/CD54 monoclonal antibody optimized for flow cytometry (FC). Detects human ICAM-1 (57kDa, SWISS P05362). Ideal for leukocyte adhesion, viral entry, and immune cell trafficking studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human ICAM-1/CD54 mAb
Remarks/Alias	ICAM1; BB2; CD54; P3.58; intercellular adhesion molecule 1
Species	Human
GeneID	3383 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-480 of human ICAM-1/CD54 (NP_000192.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P05362
Protein Weight	57kDa
Shipping	Ice bag