


## Anti-Human Ifn-Alpha/Beta R2 Monoclonal Antibody For Flow Cytometry - P48551

Item Number: CM0004502



### Introduction

-conjugated rabbit monoclonal antibody targeting human IFN-alpha/beta R2, validated for flow cytometry. Enables detection of type I interferon receptor expression in blood cells and cell lines for immunology research.

[Learn More](#)

| Parameter        | Value                                                                                                                        |
|------------------|------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human IFN-alpha/beta R2 mAb                                                                                      |
| Remarks/Alias    | IFN-R; IMD45; IFNABR; IFNARB; IFN-R-2; IFN-alpha-REC; IFNAR2                                                                 |
| Species          | Human                                                                                                                        |
| GeneID (Human)   | 3455                                                                                                                         |
| GeneID           | 3455                                                                                                                         |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 28-237 of human IFN-alpha/beta R2(NP_008896.2) |
| Source           | Rabbit                                                                                                                       |
| Category         | Monoclonal Antibodies                                                                                                        |
| Application      | FC                                                                                                                           |
| Cross Reactivity | Human                                                                                                                        |
| SWISS            | P48551                                                                                                                       |
| Protein Weight   | 58kDa/37kDa/27kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                                      |