

Anti-Human Ifn-Gamma Monoclonal Antibody For Intracellular Flow Cytometry - P01579

Item Number: CM0012118



Introduction

-conjugated rabbit monoclonal antibody targeting human IFN-gamma, validated for intracellular flow cytometry (FC). Ideal for studying type II interferon in innate and adaptive immune responses.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IFN-gamma mAb
Remarks/Alias	IFG; IFI; IMD69
Species	Human
GeneID (Human)	3458
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-166 of human IFN-gamma (NP_000610.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P01579
Protein Weight	19kDa
Shipping	Ice bag