


Anti-Human Ig Kappa Chain Monoclonal Antibody For Wb, Elisa - P01834

Item Number: CM0003413



Introduction

Recombinant rabbit monoclonal antibody against Human Ig Kappa Chain, validated for Western blotting (WB) and ELISA. Recognizes the constant region of human immunoglobulin kappa light chain, essential for antibody structure and humoral immunity studies.

[Learn More](#)

Parameter	Value
Product Name	Human Ig Kappa Chain Rabbit mAb
Remarks/Alias	Km; HCAK1; IGKCD; Human Ig Kappa Chain
Species	Human
GeneID (Human)	3514
GeneID	3514
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-107 of human Ig Kappa Chain (P01834)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P01834
Protein Weight	12kDa
Shipping	Ice bag