


Anti-Human IgA Rabbit Monoclonal Antibody For Wb And Elisa - P01876 / P01877

Item Number: CM0003240



Introduction

Rabbit monoclonal antibody against human IgA constant region (P01876) validated for Western blot and ELISA. Specific to human IgA1, no cross-reactivity with other Ig classes. Ideal for humoral immunity research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Human IgA Rabbit mAb
Remarks/Alias	IgA1; Human IgA
Species	Human
GeneID	3493/3494
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-353 of human IgA (P01876).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P01876/P01877
Protein Weight	60kDa
Shipping	Ice bag