


Anti-Human IgA1 Monoclonal Antibody () For Flow Cytometry - P01876

Item Number: CM0011010



Introduction

This -conjugated rabbit monoclonal antibody targets human IgA1 constant region, validated for flow cytometry. Ideal for detecting IgA1-expressing B cells and secreted antibodies in mucosal immunity research.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human IgA1 mAb
Remarks/Alias	IgA1
Species	Human
Gene ID (Human)	3493
Gene ID	3493
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 99-353 of human IgA1
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P01876
Protein Weight	38kDa/42KDa
Shipping	Ice bag