

Anti-Human IgD Monoclonal Antibody For Wb, Fc, Elisa - P01880

Item Number: CM0007324



Introduction

Rabbit monoclonal antibody targeting human immunoglobulin heavy constant delta (IgD). Validated for Western blot, flow cytometry, and ELISA. Suitable for B-cell biology, adaptive immunity, and cytokine induction research.

[Learn More](#)

Parameter	Value
Product Name	Human IgD Rabbit mAb
Remarks/Alias	IgD; Human IgD
Species	Human
GeneID (Human)	3495
GeneID	3495
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 99-384 of human IgD
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human
SWISS	P01880
Protein Weight	42kDa
Shipping	Ice bag