


## Anti-Human IgG (Fc) Monoclonal Antibody For IHC-P, Elisa - P01857 / P01859 / P01860 / P01861

Item Number: CM0002342



### Introduction

Mouse monoclonal antibody targeting human IgG Fc region, validated for IHC-P and ELISA. Detects human IgG heavy chain constant region (UniProt P01857, P01859, P01860, P01861); molecular weight 52kDa.

[Learn More](#)

| Parameter        | Value                                                                      |
|------------------|----------------------------------------------------------------------------|
| Product Name     | human IgG (Fc) Mouse mAb                                                   |
| Remarks/Alias    | COB1; YAP; YAP2; YAP65; YKI; YAP1; human IgG (Fc)                          |
| Species          | Human                                                                      |
| Gene ID          | 3500                                                                       |
| Immunogen        | Synthetic Peptide A synthesized peptide derived from human human IgG (Fc). |
| Source           | Mouse                                                                      |
| Category         | Monoclonal Antibodies                                                      |
| Application      | IHC-P, ELISA                                                               |
| Cross Reactivity | Human                                                                      |
| SWISS            | P01857,P01859,P01860,P01861                                                |
| Protein Weight   | 52kDa                                                                      |
| Shipping         | Ice bag                                                                    |