


Anti-Human IgG1 (Fc) Monoclonal Antibody For Wb, If-P, Elisa - P01857 / P01859 / P01860 / P01861

Item Number: CM0003723



Introduction

Rabbit monoclonal antibody targeting human IgG1 Fc region. Recombinant immunogen, recommended for WB, IF-P, and ELISA. Cross-reacts with human. Facilitates study of immunoglobulin constant region functions in adaptive immunity.

[Learn More](#)

Parameter	Value
Product Name	Human IgG1 (Fc) Rabbit mAb
Species	Human
GeneID (Human)	3500
GeneID	3500
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-330 of human IgG1 (Fc) (P01857)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human
SWISS	P01857, P01859, P01860, P01861
Protein Weight	44kDa
Shipping	Ice bag