


Anti-Human IgG3 Polyclonal Antibody For Wb And Elisa - P01860

Item Number: CM0020371



Introduction

Rabbit polyclonal antibody targeting the human IgG3 heavy chain constant region (aa 80-180), validated for Western blot and ELISA. Ideal for investigating humoral immunity, antibody effector mechanisms, and B-cell biology.

[Learn More](#)

Parameter	Value
Product Name	Human IgG3 Rabbit pAb
Remarks/Alias	IgG3; Human IgG3
Species	Human
Gene ID (Human)	3502
Gene ID	3502
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-180 of human IgG3. (P01860)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P01860
Protein Weight	41kDa
Shipping	Ice bag