


Anti-Human IgG3 Rabbit Monoclonal Antibody For Wb, Elisa - P01860

Item Number: CM0007519



Introduction

Human IgG3 rabbit monoclonal antibody suitable for WB and ELISA. Targets the constant region (aa 99-167) of human IgG3 (UniProt P01860) for humoral immunity research and antibody detection.

[Learn More](#)

Parameter	Value
Product Name	Human IgG3 Rabbit mAb
Remarks/Alias	IgG3; Human IgG3
Species	Human
Gene ID (human)	3502
Gene ID	3502
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P01860
Protein Weight	41KDa/49KDa
Shipping	Ice bag