


Anti-Human IgG4/Immunoglobulin Heavy Constant Gamma 4 Rabbit Polyclonal Antibody For Wb And Elisa - P01861

Item Number: CM0018876



Introduction

Rabbit polyclonal antibody targeting human IgG4 constant region (IGHG4). Validated for Western blot and ELISA. Recognizes Ig gamma-4 chain C region for humoral immunity and immune tolerance research.

[Learn More](#)

Parameter	Value
Product Name	Human IgG4/Immunoglobulin heavy constant gamma 4 Rabbit pAb
Remarks/Alias	IGHG4; Ig gamma-4 chain C region; Human IgG4/Immunoglobulin heavy constant gamma 4
Species	Human
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human IgG4/Immunoglobulin heavy constant gamma 4 (P01861)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P01861
Protein Weight	44kDa
Shipping	Ice bag