


Anti-Human IgM (Fc) Rabbit Monoclonal Antibody For Wb, If/Icc, Fc, Elisa - P01871

Item Number: CM0006386



Introduction

Recombinant rabbit monoclonal antibody targeting the constant region of human IgM heavy chain (P01871). Suitable for ELISA, WB, IF/ICC, and FC applications. Ideal for humoral immunity, BCR signaling, and IgM effector function studies.

[Learn More](#)

Parameter	Value
Product Name	Human IgM (FC) Rabbit mAb
Remarks/Alias	MU; VH; AGM1; Human IgM
Species	Human
GeneID (Human)	3507
GeneID	3507
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 218-453 of Human IgM (P01871).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, FC, ELISA
Cross Reactivity	Human
SWISS	P01871
Protein Weight	49kDa
Shipping	Ice bag