


Anti-Human Igm Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P01871

Item Number: CM0004250



Introduction

Rabbit monoclonal antibody specific for Human IgM, suitable for WB, IHC-P, IF/ICC, ELISA. Recognizes the constant region of human IgM heavy chain (~49kDa). Ideal for B cell and humoral immunity research.

[Learn More](#)

Parameter	Value
Product Name	Human IgM Rabbit mAb
Remarks/Alias	MU; VH; AGM1; Human IgM
Species	Human
GeneID (human)	3507
GeneID	3507
Immunogen	`Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P01871
Protein Weight	49kDa
Shipping	Ice bag