


Anti-Human IL-17A Monoclonal Antibody For Intracellular Flow Cytometry - Q16552

Item Number: CM0005149



Introduction

-conjugated anti-human IL-17A monoclonal antibody validated for intracellular flow cytometry. Detects the pro-inflammatory cytokine IL-17A (SWISS: Q16552), expressed in Th17 cells. Suitable for studying host defense, inflammation, and autoimmune mechanisms.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IL-17A mAb
Remarks/Alias	IL17; CTLA8; IL-17; ILA17; CTLA-8; IL-17A
Species	Human
GeneID (human)	3605
GeneID	3605
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-155 of human IL-17A (NP_002181.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	Intracellular Flow Cytometry (FC)
Cross Reactivity	Human
SWISS	Q16552
Protein Weight	18kDa
Shipping	Ice bag