

Anti-Human IL-17Rb Monoclonal Antibody () For Flow Cytometry - Q9Nrm6

Item Number: CM0012655



Introduction

-conjugated rabbit monoclonal antibody against human IL-17RB (Q9NRM6) for flow cytometry. Ideal for profiling IL-17B/E receptor expression in endocrine tissues and hematopoietic research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IL-17RB mAb
Remarks/Alias	IL17RB; CRL4; EVI27; IL17BR; IL17RH1; interleukin-17 receptor B
Species	Human
Gene ID	55540
Immunogen	Recombinant fusion protein containing amino acids 18-289 of human IL-17RB (NP_061195.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9NRM6
Protein Weight	31kDa/55kDa
Shipping	Ice bag