

Anti-Human IL-1B Monoclonal Antibody () For Intracellular Flow Cytometry - P01584

Item Number: CM0011201



Introduction

-conjugated rabbit anti-human IL-1 β monoclonal antibody for intracellular flow cytometry. Detects human IL-1 β , a key pro-inflammatory cytokine. Ideal for inflammation and immune response research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IL-1 β mAb
Remarks/Alias	IL1B; IL-1; IL1-BETA; IL1F2; interleukin-1 beta
Species	Human
GeneID (Human)	3553
GeneID	3553
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 117-269 of human IL1 β (NP_000567.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P01584
Protein Weight	30kDa
Shipping	Ice bag