

Anti-Human IL-4 Monoclonal Antibody () For Intracellular Flow Cytometry - P05112

Item Number: CM0009701



Introduction

High-specificity -conjugated rabbit anti-human IL-4 monoclonal antibody, validated for intracellular flow cytometry. Ideal for Th2 immune profiling, cytokine detection in allergy, asthma, and immuno-oncology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IL-4 mAb
Remarks/Alias	BSF1; IL-4; BCGF1; BSF-1; BCGF-1
Species	Human
GeneID (Human)	3565
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-153 of human IL-4 (NP_000580.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P05112
Protein Weight	17 kDa
Shipping	Ice bag