


Anti-Human IL-5 Monoclonal Antibody For Flow Cytometry - P05113

Item Number: CM0009189



Introduction

-conjugated monoclonal antibody targeting human IL-5 (Interleukin-5). Suitable for intracellular flow cytometry (FC). Research use for eosinophil differentiation, B-cell activation, and inflammatory diseases. SWISS: P05113.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human IL-5 mAb
Remarks/Alias	EDF; IL-5; TRF; interleukin-5
Species	Human
GeneID (Human)	3567
GeneID	3567
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-134 of human IL-5 (NP_000870.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P05113
Protein Weight	15kDa
Shipping	Ice bag