


Anti-Human Ilt5/Cd85A Monoclonal Antibody For Flow Cytometry - O75022

Item Number: CM0009336



Introduction

Rabbit anti-Human ILT5/CD85a monoclonal antibody validated for flow cytometry (FC). Recognizes human ILT5/CD85a, an inhibitory immune receptor on monocytes and B-cells. Ideal for immunophenotyping and functional studies in immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human ILT5/CD85a mAb
Remarks/Alias	LILRB3; CD85A; HL9; ILT-5; ILT5; LILRA6; LIR-3; LIR3; PIR-B; PIRB; leukocyte immunoglobulin like receptor B3
Species	Human
GeneID (Human)	11025
GeneID	11025
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-418 of human ILT5/CD85a (NP_006855.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O75022
Protein Weight	69kDa/71kDa
Shipping	Ice bag