


Anti-Human Integrin B5/Itgb5 Rabbit Monoclonal Antibody () For Flow Cytometry - P18084

Item Number: CM0008235



Introduction

This -conjugated rabbit monoclonal antibody targets human Integrin $\beta 5$ /ITGB5, validated for flow cytometry. Ideal for studying fibronectin adhesion, integrin signaling, and adenovirus receptor function.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Integrin $\beta 5$ /ITGB5 mAb
Remarks/Alias	ITGB5; integrin beta-5
Species	Human
GeneID (Human)	3693
GeneID	3693
Immunogen	Recombinant protein corresponding to a sequence within amino acids 24-490 of human Integrin $\beta 5$ /ITGB5 (NP_002204.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P18084
Protein Weight	88kDa
Shipping	Ice bag